

Respirator

FFP2



Features

Material	PP
Inner Layer	Hydrophobic / Hydrophilic
Ply	3
Tie Type	Straps
Valve	With / Without
Nose Bridge	Enclosed / Exposed
Sterilization	Non-sterile

* Tested by Nelson Labs

Standard

▶ EN 149

Colors



White



Navy Blue

▶ Customized print is also available

Size

160 x 110 mm

Packaging

50 pcs / box

10 boxes / case

Ordering Code

GL001 (With Valve)

GL001A (Without Valve)