

# Respirator

KN95

NIOSH

FFP2



## Features

Material	PP
Inner Layer	Hydrophobic / Hydrophilic
Ply	3 / 4
Tie Type	Earloops / Straps
Valve	With / Without
Nose Bridge	Enclosed / Exposed
Sterilization	Sterile / Non-sterile

\* Tested by Nelson Labs

## Standards

- ▶ GB 2626
- ▶ NIOSH
- ▶ EN 149

## Colors



White



Navy Blue

- ▶ Customized print is also available

## Size

160 x 110 mm

## Packaging

## Ordering Code

GB 2626

40 pcs / box  
30 boxes / case

D918

NIOSH

50 pcs / box  
10 boxes / case

901AV

EN 149

50 pcs / box  
10 boxes / case

GL001 (With Valve)  
GL001A (Without Valve)