



Item	Type I+ Procedure Mask , Green
Standard	(Type I+) BFE > 95%
Sterile	Non-sterile
Description	Airborne protective equipment
Position of Loop	Outside
Outer layer (face out)	20gsm PP
Middle Layer(Filter)	22gsm meltblown
Inner layer (face in)	30gsm ES
Ear-loop	Round
Nose Bridge	Metal stripe coated by PE
Color	Green
Dimensions	Adult Size 175±2mm x 95±2mm
Packing	Bulk Pack: 50pcs/box 20 boxes / case Box dimension : 18x10x8cm Case dimension: 37x42x21cm
Antifog	/

