



Item	Type II Antifog Procedure Mask , Blue
Standard	(Type II) BFE > 98%
Sterile	Non-sterile
Description	Airborne protective equipment
Position of Loop	Outside
Outer layer (face out)	16gsm PP
Middle Layer (Filter)	25gsm meltblown
Inner layer (face in)	30gsm ES
Ear-loop	Flat
Nose Bridge	Metal stripe coated by PE
Color	Blue
Dimensions	Adult Size 175±2mm x 95±2mm
Packing	Bulk Pack: 50pcs/box , 20 boxes / case Individual Pack : 1pcs/Bag ,50 Bags/ Box ,20 boxes / case Box dimension : 18x10x8cm / 20x10x13cm Case dimension: 37x42x21cm
Antifog	Min

