



Item	Type II Infant Mask , Blue
Standard	( Type II ) BFE > 98%
Sterile	Non-sterile
Description	Airborne protective equipment
Position of Loop	Outside
Outer layer (face out)	20gsm PP
Middle Layer( Filter)	25gsm meltblown
Inner layer (face in)	30gsm ES
Ear-loop	Flat
Nose Bridge	Metal stripe coated by PE
Color	Blue
Dimensions	Size $120 \pm 2\text{mm} \times 85 \pm 2\text{mm}$
Packing	Bulk Pack: 50pcs/box , 20 boxes / case Box dimension : 8.7x12.4x10cm Case dimension: 27x47.5x22cm

