



| | |
|------------------------|--|
| Item | Type II+ Procedure Mask , Blue |
| Standard | (Type II+) BFE > 98% |
| Sterile | Non-sterile |
| Description | Airborne protective equipment |
| Position of Loop | Outside |
| Outer layer (face out) | 20gsm PP |
| Middle Layer (Filter) | 25gsm meltblown |
| Inner layer (face in) | 30gsm ES |
| Ear-loop | Flat |
| Nose Bridge | Metal stripe coated by PE |
| Color | Blue |
| Dimensions | Adult Size 175±2mm x 95±2mm |
| Packing | Bulk Pack: 50pcs/box 20 boxes / case Box dimension : 18x10x8cm Case dimension: 37x42x21cm |
| Antifog | / |

