



Item	Type IIR 160 Antifog Procedure Mask , Blue
Standard	( Type IIR ) BFE > 98%
Sterile	Non-sterile
Description	Airborne protective equipment
Position of Loop	Outside
Outer layer (face out)	25gsm PP
Middle Layer( Filter)	25gsm meltblown
Inner layer (face in)	40gsm ES
Ear-loop	Flat
Nose Bridge	Metal stripe coated by PE
Color	Blue
Dimensions	Adult Size 175±2mm x 95±2mm
Packing	Bulk Pack: 50pcs/box 20 boxes / case Box dimension : 18x10x8cm Case dimension: 37x42x21cm
Antifog	Max

