



Item	Type IIR 160 SURGICAL MASK Antifog, Tie-On
Standard	(Type IIR) BFE > 98%
Sterile	Non-sterile
Description	Airborne protective equipment
Position of Loop	Outside
Outer layer (face out)	25gsm PP
Middle Layer (Filter)	25gsm meltblown
Inner layer (face in)	40gsm ES
Ear-loop	Tie
Nose Bridge	Metal stripe coated by PE
Color	Blue
Dimensions	Adult Size 175±2mm x 95±2mm
Packing	Bulk Pack: 50pcs/box , 40 boxes / case Box dimension : 18.5x10.5x1.5cm Case dimension: 56.5x39.5x49cm
Antifog	Max

